

# STERI-PEROX 3%

SPER-16Z-3

## DESCRIPTION

STERI-PEROX 3% was developed for use in cleaning rotation cycles to maintain a clean environment and reduce the regulatory concern for VOCs. STERI-PEROX 3% is an effective one step ready to use solution. The product is packaged to comply with the requirements for transportation into a classified area. 16oz.

## FEATURES & BENEFITS

- » Each sterile container is double bagged packaged
- » Quadruple bagged in the ABCD Cleanroom Introduction System®
- » Individually labeled with lot number and expiration
- » Delivered with lot specific Certificate of Analysis and Certificate of Sterility
- » Specifically formulated as a sterile cleanroom pharmaceutical formula
- » Compatible with most surfaces
- » 16 oz containers comes with sterile spray attachment
- » Low remaining residue

## QUALITY & MANUFACTURING

- » Formulated with USP Water for Injection
- » Filled in an ISO 5 (Grade A/B, Former Class 100)
- » Filtered at 0.2 microns
- » Raw components are air washed with 0.2 micron filtered air
- » Aseptically filled into sterile components via gamma irradiation
- » Lot sterility tested according to current USP compendium
- » Completely traceable from start to finish
- » Completely validated for sterility and shelf life

## PACKAGING

12/case

